M-30-15-0685(04)

User's Manual
ENGLISH

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Promedon

People + Innovation

USER'S MANUAL

DESCRIPTION

Safyre is a kit for the treatment of stress urinary incontinence grades II and III and it contains:

- 1 **Safyre** Sling (Ref: S-85-AF), made out of highly biocompatible synth materials.
- 2 Deschamps needles (Ref: DPN-HA), This single-use surgical instrument is intended to be used in combination with the Sling for its implantation.
- 2 pairs of radiopaque Washers (AR/TP).

All components of **Safyre** are sold sterile and ready to use.

Safyre Sling is a permanent implant that consists of two columns formed by multiple cone-shaped subunits made of polydimethylsiloxane and a central polypropylene mesh that surrounds the urethra under the mid urethra, producing a coaptation of the former during strains. The columns formed by cones autofix to the obturator muscle anchoring the sling. retropubic spaces anchoring the sling.

CONTRAINDICATIONS

- **Safyre** cannot be prescribed if there is some kind of infection, specially genital or related to the urinary tract.
- **Safyre** must not be used with patients under anticoagulation therapy, during pregnancy or with an existing urinary tract infection.
- Do not use **Safyre** the patient has cystocele grade III.
- **Safyre** should be used with precaution in diabetic women.

WARNING

The prosthesis must not be manipulated with pointed, serrated or sharp elements since any hollow, damage, hole or tear can be the cause of subsequent complications. Fluffs, fingerprints, talcum powder or any contaminated surface can be the cause of reactions to foreign bodies. Maximum precautions should be taken to avoid contamination.

Safyre sling must be implanted without tension, that is, the central mesh should rest on the urethra without any tension. An implantation with tension leads to problems of urethral erosion and urinary retention. Punctures or injuries of blood vessels, bowels or nerves could occur when passing the needle and a surgical repair may be needed.

As with all foreign bodies, the polypropylene mesh and silicone columns could potentiate an existing infection. The wrapping around the sling is intended to minimize the risk of contamination during the surgical procedure. The components of **Safyre** have been designed to be used only ONCE. Therefore DO NOT RE-USE OR RE-STERILIZE.

INFORMATION FOR THE PATIENT

The surgeon has the responsibility of informing the patient or her representative about the possible complications related to the implantation of **Safyre** Sling before the surgery. **Promedon** and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and the possible risks related to the implantation and use of the Sling. The patient should be advised that the implantation of the sling might loose effectiveness with future pregnancies, and that she could become incontinent again. It is convenient that the patient avoids weight lifting or intense exercises that involve effort (cycling, aerobics, etc.) for the first three to four weeks after the surgery, and sexual relationships for at least one month.

It is essential that the patient contacts the surgeon immediately in case of:

- Dysuria (pain or difficulty during urination).
- Vaginal pain.
- Fever.
- Presence of serous, bleeding or purulent secretion.
- Bleeding or other problems.

PRECAUTIONS

A careful patient selection as well as a complete diagnosis previous to surgery is essential.

The sling should be carefully manipulated avoiding the use of pointed, serrated

or sharp elements.

Fluffs, fingerprints, talcum powder, bacteria or other elements inclined to contaminate the surface of sling may be the cause of infections or reactions to foreign bodies. Maximum precautions should be taken to avoid contamination. That is why it is suggested to keep the sling in contact with antibiotics during surgery.

The use of a sterile paste composed of 10 grams of Xylocaine Gel and two ampoules (160 mg/ampoule) of Gentamycine is specially recommended. This mixture forms a gelatinous paste that adheres to the sling, and when implanting it part of the mixture enters the body.

The surgical procedure should be performed with care to avoid large vessels, nerves, and bowels. Attention to local anatomy and proper passage of needle will minimize risks.

Manipulation and Storage Precautions:  is provided sterile and pyrogen-free. The packaging consists of a cardboard box that contains two envelopes: one with a  Sling and the other with the surgical instrument. IF AN ENVELOPE IS DAMAGED, DO NOT IMPLANT THE SLING.

Operating Room Conditions:

In general, all operating rooms fulfill the necessary conditions to implant the sling. Nevertheless, it is advisable that the following parameters are specially checked:

- Asepsis of the Operating Room.
- Adequate Training of the Operating Room Assistant Personnel.

STORAGE

 should be stored under the following conditions:

- **TEMPERATURE:** ROOM TEMPERATURE

DO NOT USE AFTER THE EXPIRATION DATE SHOWN IN THE PACKAGING.

COMPLICATIONS

The possible complications associated with the use of sling must be discussed with the patient before surgery.

The use of this sling may result in some complications associated with the medication and the methods used in the surgical procedure as well as the reaction of the patient or the degree of intolerance to any foreign body implanted in the body. Some complications may require the sling removal. Infections that do not respond to an antibiotic therapy require the prosthesis removal. Some patients may experience vaginal pain during the initial post-surgery period. A therapy with ANALGESIC and ANTI-INFLAMMATORY drugs could be enough to relieve the pain.

Other complications reported with this or other slings include the following:

- Wound infection including a secondary necrosis.
- A serious infection cannot be treated with the prosthesis implanted.
- Urethral or vaginal erosion.
- Vaginal pain.
- Serous or blood purulent flow.
- Vaginal inflammation.
- Injuries of vessels and nerves.
- Presence of vaginal fistula.
- Vesical instability.
- Vaginal Dehiscence.
- Urinary obstruction.

Post-surgery formation of a capsule of fibrous tissue around the sling is a normal physiologic answer to the implantation of a foreign body.

In case of vaginal exteriorization of the sling, often due to infection, it is necessary to remove the implant. This can be achieved by simply cutting the polypropylene mesh and pulling both columns out.

All surgeons must notify **Promedon** about any complication developed by the use of .

SURGICAL PROCEDURE

Prepare the patient for the surgery as usual, and introduce a Foley catheter in the urethra. During the surgery, we recommend keeping the sling in its inner wrapping (sheath).

 Sling is implanted in a transobturator approach, under regional or local anesthesia, following the surgery technique currently accepted and used by the present art. The routine use of prophylactic IV antibiotics should be considered.

A description of the transobturator approach is summarized in the following steps:

1. Medium Colpotomy: Make a sagittal incision about 1.5cm long starting approximately 1cm from the outer urethral meatus.

2. Paraurethral Dissection: From the incision, exteriorize the vaginal wall with a pair of scissors. Minimal vaginal dissection is needed to achieve a tunnel that allows the passage of the transobturator needle.

3. Sling Placement: Bilaterally make a punctiform skin incision in the genitofemoral folds at the level of the clitoris. The incision should be made at the same horizontal level of the clitoris.

Needles are introduced parallel to the ischiopubic ramus and exteriorized through the vaginal incision, simply by fists rotation, guided by the surgeon’s forefinger.

Thread the tip of the needle through the perforation in the column of the sling and transfer the column through the tunnel previously created by the needle. If tear at the perforation of the sling column while it is being transferred through the patient’s tissues occurred, tie a suture at the end of the column. Threading the suture in the hole of the needle shall allow completing the transference of the sling safely. Repeat this step over the other obturator hole.

4. Tension-Free Fixation: Place a pair of Metzenbaum scissors between the mesh and the urethra to facilitate the regulation of tension and to avoid the mesh folding. Gently pull up the columns to bring the mesh in contact with the urethra. The right tension can be determined with an intraoperative cough stress test verifying continence. Pull the columns up to increase the tension on the urethral support, or release them to loosen it. (Do not loosen the sling by the polypropylene mesh). Finally, cut the exceeding parts of the columns and suture the skin incisions.

Post-operative cares and therapy are left to the Surgeon’s criteria.

Suggestion 1: We suggest cutting the column 7 cones above the skin and placing the column excess tunneled in the subcutaneous cellular tissue. This practice is particularly beneficial, providing for better locating the columns when making post-surgical readjustments.

Suggestion 2 : The use of washers is optional and depends on surgeon’s criteria. It is recommended for those cases in which the sling must be implanted with some tension, as in the case of rigid or poor urethras.

Washer placement: Prior to placing the washers, tissues should be extended with the use of forceps in order to create space for placing them. The washer is threaded into mosquito forceps which will then be used to hold the end of the column in order to transfer the washer from the forceps to the column. Then, the washer is moved down the column to the obturator wall. Finally, cut the remainder of the column verifying that there are at least 3 cones above the washer.

Washers are radiopaque, which makes it easier to locate the end of the column in the case of postoperative readjustment.

SYMBOLS USED IN THE LABELS

CATALOGUE NUMBER

BATCH NUMBER

SEE INSTRUCTIONS FOR USE

EXPIRATION DATE

DO NOT REUSE

DATE OF MANUFACTURE

STERILE: STERILIZATION METHOD
ETHYLENE OXIDE

MANUFACTURER

AUTHORIZED REPRESENTATIVE
IN THE EUROPEAN COMMUNITY

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